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Frozen slides experimental report

1. Experimental equipments and reagents

1.1 Experimental equipments

Name	Manufacturer	Model
Frozen section machine	Shenyang Deyu Electronic Instrument Co., Ltd	SYD-K2040
Adhesion microscope slide	Shaanxi Yike Biotechnology Co., Ltd	YK116
Slicer	Shanghai Leica Instrument Co., Ltd	LEICA 819

1.2 Main experimental reagents

Reagent	Manufacturer	Article number
OCT embedding agent	SAKURA (樱花)	4583
Polyformaldehyde fixsolution	Shaanxi Yike Biotechnology Co., Ltd	YK2200
Sucrose	Sinopharm Group Chemical Reagent Co. LTD	57-50-1

2. Frozen slides making

2.1 Tissue fixation: Fresh tissues that need to be fixed were fixed in the fixative for more than 24 hours. However, the fresh and directly sliced tissues stored in low temperature transportation skip the steps of fixation and dehydration, and are embedded and sliced by Oct.

2.2 Dehydration: The tissue was removed from the fixed solution and placed in 15% sucrose solution for 4 °C refrigerator dehydration, then transferred into 30% sucrose solution for 4 °C refrigerator dehydration.

2.3 OCT embedding: Take out the dehydrated tissue, drain the surface water with filter paper, smooth the target tissue with scalpel, put the section face up on the specimen chuck, drop OCT embedding agent around the tissue, put the specimen chuck on the quick freezing table of the frozen section machine for quick freezing and embedding. When OCT becomes

white and hard, it can be sectioned.

2.4 Section: The specimen chuck was fixed on the slicer, the tissue surface was rough cut, and then sliced. The slice thickness was 8-10 μ M. the clean slide was placed on the top of the cut tissue slide, and the tissue was pasted on the slide. After the film is labeled, it should be stored at - 20 °C for future use.