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---

Fluorescent-labelled Secondary antibody	Xi'an ambition Biotechnology Co., LTD
Spontaneous fluorescence quenching reagent	BOSTER
DAPI	MCE

---

## 2 Procedure

**2.1 Deparaffinize and rehydrate:** incubate sections in 3 changes of Biodewax and Clear Solution, 10 min each. Dehydrate in 3 changes of pure ethanol for 5 min. Wash in distilled water.

**2.2 Antigen retrieval:** The tissue sections were placed in a pressure cooker with citric acid (PH 6.0) antigen repair solution for antigen repair. Spray time for 2.5min. After natural cooling, the slides were placed in PBS (PH 7.4) and washed for 3 times for 5min each time.

**2.3 Circle and Serum blocking:** eliminate obvious liquid, mark the objective tissue with liquid blocker pen. Add 5% BSA to cover the marked tissue to block non-specific binding for 30 min. Cover objective area with 10% donkey serum (for the case of primary antibody originated from goat) or 5% BSA (for the case of primary antibody originated from others).

**2.4 Primary antibody:** throw away the blocking solution slightly. Incubate slides with primary antibody (diluted with PBS appropriately) overnight at 4 °C, placed in a wet box containing a little water.

**2.5 Secondary antibody:** wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then throw away liquid slightly. Cover objective tissue with secondary antibody (appropriately respond to primary antibody in species), incubate at room temperature for 50 min in dark condition.

**2.6 DAPI counterstain in nucleus:** wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then incubate with DAPI solution at room temperature for 10 min, kept in dark place.

**2.7 Spontaneous fluorescence quenching :** wash three times with PBS (pH 7.4) in a Rocker device, 5min each. Add spontaneous fluorescence quenching reagent to incubate for 5 min. Wash in running tap water for 10 min.

**2.8 Mount:** Throw away liquid slightly, then cover slip with anti-fade mounting medium.

**2.9 Microscopy detection and collect images by Fluorescent Microscopy.** DAPI glows blue by UV excitation wavelength 330-380nm and emission wavelength 420 nm; FITC glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm.

## 3 Results

Nucleus is blue by labeling with DAPI. Positive cells are green or red according to the fluorescent labels used.