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## Paraffin sides experiment report

---

### 1. Experimental equipment and reagents

#### 1.1 Experimental equipment

Name	Manufacturer	Model
Dehydrator	Wuhan Junjie Electronics Co., Ltd	JJ-12J
Embedding machine	Wuhan Junjie Electronics Co., Ltd	JB-P5
Pathology slicer	Shanghai Leica Instrument Co., Ltd	RM2016
Frozen platform	Wuhan Junjie Electronics Co., Ltd	JB-L5
Tissue spreader	Zhejiang Jinhua Kdy Equipment Co., Ltd	KD-P
Oven	Shanghai Huitai Instrument Manufacturing Co., Ltd	DHG-9140A
Slides	Shaanxi Yike Biotechnology Service Co., Ltd	YK116

#### 1.2 Main experimental reagents

Reagent	Manufacturer	Article number
Xylene	Sinopharm Group Chemical Reagent Co. Ltd	10023418
Ethanol	Sinopharm Group Chemical Reagent Co. Ltd	100092683

### 2. Paraffin section making

**2.1 draw materials:** The fresh tissue was fixed with fixed liquid for more than 24 hours. Remove the tissue from the fixed liquid and trim the target tissue with a scalpel in the ventilation cupboard, and put the trimmed tissue and the label in the dehydration box.

**2.2 Dehydration and wax leaching:** put the dehydration box into the dehydrator in order to dehydrate with gradient alcohol. 75% alcohol 2 hours, 85% alcohol 1 hours, 90% alcohol 1 hours, 95% alcohol 1 hour, anhydrous ethanol II h, anhydrous ethanol II 1h, xylene I 1 h, xylene II 1 h, 56-58°C melting paraffin I 1h, 56-58°C melting paraffin II 1h, 56-58°C melting paraffin III 1 hour.

**2.3 embedding:** The wax-soaked tissue is embedded in the embedding machine. First, put the melted wax into the embedding frame, and before the wax solidifies, remove the tissue from the

---

dewatering box and put it into the embedding frame according to the requirements of the embedding surface and affix the corresponding label. Cool at -20 °C freezing table, and after the wax is solidified, the wax block is removed from the embedded frame and repaired.

**2.4 cut into slices:** Place the trimmed wax block cool at -20 °C freezing table, slice the modified tissue chip wax block on the paraffin slicer, the slice thickness is 4 μm. The tissue is flattened when the slice floats on the 40 °C warm water of the spreading machine, and the tissue is picked up by the glass slides and baked in the oven at 60 °C. After the water-baked dried wax is melted, it is taken out and stored at room temperature.